



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vapour Pressure Of Water Marcet Boiler

Product Code :
LBNY-0005-103000125



Description :

Measuring a vapour pressure curve for saturated vapour
Bourdon tube pressure gauge to indicate pressure
Digital temperature display
Boiler with insulating jacket
Temperature limiter and safety valve protect against overpressure in the system.

Technical Specification :

Technical Data:
Bourdon tube pressure gauge: -1...24bar
Heater: 2kW
Safety valve: 20bar
Temperature limiter: 200C
Boiler, stainless steel: 2L
Measuring ranges
Pressure: 0...20bar
Temperature: 0...250C.

Laboratory Instrument India

