



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heat Exchanger Supply Unit

Product Code :
LBNY-0005-103000126



Description :

Supply unit for heat exchangers
Digital displays for 6 temperature and 2 flow rate sensors
Flow adjustable using valves
Water connections with quick-release couplings
Hot water circuit with tank, heater, temperature controller, pump and protection against lack of water
Temperature controller controls the temperature of hot water
Stirring machine connection with speed adjustment
Cold water circuit from laboratory mains or water chiller.

Technical Specification :

Technical Data:

Heater
Power output: 3kW
Hot water tank: approx. 10L
Thermostat: 0...70C
Pump
Max. flow rate: 600L/h
Power consumption: 120W
Max. head: 30m
Measuring ranges
Flow rate: 2x 20...250L/h
Temperature: 6x 0...100C.



Laboratory Instrument India