



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Heat Exchanger Supply Unit

**Product Code :**  
LBNY-0005-103000126



#### **Description :**

Supply unit for heat exchangers  
Digital displays for 6 temperature and 2 flow rate sensors  
Flow adjustable using valves  
Water connections with quick-release couplings  
Hot water circuit with tank, heater, temperature controller, pump and protection against lack of water  
Temperature controller controls the temperature of hot water  
Stirring machine connection with speed adjustment  
Cold water circuit from laboratory mains or water chiller.

#### **Technical Specification :**

##### **Technical Data:**

Heater  
Power output: 3kW  
Hot water tank: approx. 10L  
Thermostat: 0...70C  
Pump  
Max. flow rate: 600L/h  
Power consumption: 120W  
Max. head: 30m  
Measuring ranges  
Flow rate: 2x 20...250L/h  
Temperature: 6x 0...100C.



**Laboratory Instrument India**