




Email : sales@laboratoryinstrumentindia.com

Product Name :
Two Stroke Petrol Engine

Product Code :
LBNY-0005-103000128



Description :

Engine completely equipped with fuel line, throttle cable and exhaust gas temperature sensor
Air-cooled single-cylinder two-stroke petrol engine for installation in test stand
Fuel hose with self-sealing quick-release coupling
Force transmission to brake via pulley, gear transmission 2:1
Engine mounted on a base plate with vibration dampers.

Technical Specification :

Technical Data:

Air-cooled single-cylinder petrol engine
Bore: 42,5mm
V-belt: diameter=63mm
Stroke: 32mm
Power output: 1,32kW at 7000min-1.



Laboratory Instrument India

