



Email : sales@laboratoryinstrumentindia.com

Product Name :
Ground Water Flow

Product Code :
LBNY-0005-103000141



Description :

Investigation of groundwater flows
19 measuring connections with filters to detect the groundwater levels, arranged orthogonal to the tank bottom
Stainless steel tank as experimental section to be filled with coarse sand
2 different models for excavation pits
Water drain via 2 wells with open-seam tubes in the experimental section
Water supply via 2 open-seam tubes
Water feeds and discharges can be adjusted separately via valves
Groundwater levels displayed on the 19 tube manometers
1 model for structure with waterproof bottom.

Technical Specification :

Technical Data:

Tube manometers
Number of tubes: 19
Measuring range: 300mmWC
Tank
Material: stainless steel
Content, LxWxH: 1000x615x350mm
19 measuring connections on the bottom of the tank
Plastic models
Excavation pit, LxWxH: 256x464x150mm
Excavation pit, LxWxH: 610x464x150mm

Structure with waterproof bottom,
DxH: 180x150mm, inner tube DxH: 40x330m.



Laboratory Instrument India