



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wave Generator

Product Code :
LBNY-0005-103000143



Description :

Generation of surface waves via paddle swinging back and forth
Motor with variable speed via frequency converter
Wave generator for the experimental flume
Drive of the paddle via crank mechanism and motor with gear
Weir at outlet element of used as paddle
Adjustable crank mechanism stroke.

Technical Specification :

Technical Data:

Motor
Stroke: 40...200mm
Power output: 0,55kW
Output speed: 0...110min-1
Speed: 1375min-1.



Laboratory Instrument India

