



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cake And Depth Filtration

Product Code :
LBNY-0005-103000145



Description :

Pressure loss measurement with twin tube manometers
Height-adjustable filler hopper made of glass
Flow meter with needle valve for adjustment
Fundamentals of cake and depth filtration
Filter element with sintered filter medium on its bottom to capture the particles

Technical Specification :

Technical Data:

Filter element
Filter chamber height: 85mm
Cross-sectional area: approx. 11cm
Inside diameter: approx. 37mm
Tube material: glass
Filter medium, sintered filter
Material: sintered metal
Pore size: 100m
Thickness: 2mm
Measuring ranges:
Pressure: 2x 0...500mmWC
Flow rate: 40...360mL/min
Temperature: -10...100C
Measuring cups

1x 100mL, scale division: 2mL
1x 1000mL, scale division: 10mL.



Laboratory Instrument India