



Email : sales@laboratoryinstrumentindia.com

Product Name :
Full Scale Sewerage System

Product Code :
LBNY-0005-103000148



Description :

6 tube manometers to indicate the pressure curve in the downcomer
Measurement of flow rate via rotameter
Experimental plant for demonstration of wastewater technology
Transparent glass pipes and tanks
Contains downpipe, collection pipe, ventilation pipe and bypasses
Ten cisterns with remotely-operated solenoid valves
Two toilets with cistern or pressure flush.

Technical Specification :

Technical Data:

Pump
Max. head: 42m
Max. flow rate: 70L/min
Power consumption: 550W
Collection tank
Volume approx. 300L
Transparent cisterns
4x 20L
6x 10L
1x 9L
Flush for toilet max. 9L
Measuring ranges:

Tube manometer, 6-fold 1500mm
Flow rate: 0,4...4L/h.



Laboratory Instrument India