



Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard Chemicals Pump

Product Code :
LBNY-0005-103000152



Description :

Process configuration permits easy exchange of wearing parts
Pump hydraulic design according
Centrifugal pump as accessory for installation
Pump technical requirements according
Drive and water supply provided.

Technical Specification :

Technical Data:

Centrifugal pump (at nominal speed: 1450min-1)
Max. head: 9,5m
Max. flow rate: 9,5m³/h
Power consumption: 0,5kW
Materials
Shaft: stainless steel
Housing, impeller: grey cast iron.



Laboratory Instrument India

