



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Standard Chemicals Pump With Magnetic Clutch

**Product Code :**  
LBNY-0005-103000153



**Description :**

Permanent-magnetic synchronous drive inside pump  
Single-stage centrifugal pump with magnetic clutch as accessory for installation  
Drive and water supply provided  
Pump technical requirements.

**Technical Specification :**

**Technical Data:**

Pump (at nominal speed: 2900min-1)  
Power consumption: 3,7kW  
Max. head: 39m  
Max. flow rate: 12m<sup>3</sup>/h  
Dimensions And Weight:  
LxWxH: 625x240x300mm  
Weight: 60kg.



**Laboratory Instrument India**

