



Email : sales@laboratoryinstrumentindia.com

Product Name :
Standard Chemicals Pump With Magnetic Clutch

Product Code :
LBNY-0005-103000153



Description :

Permanent-magnetic synchronous drive inside pump
Single-stage centrifugal pump with magnetic clutch as accessory for installation
Drive and water supply provided
Pump technical requirements.

Technical Specification :

Technical Data:

Pump (at nominal speed: 2900min-1)
Power consumption: 3,7kW
Max. head: 39m
Max. flow rate: 12m³/h
Dimensions And Weight:
LxWxH: 625x240x300mm
Weight: 60kg.



Laboratory Instrument India

