



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Backflow Preventer

Product Code :
LBNY-0005-103000158



Description :

Functionality of moving parts preserved
Unit designed for desktop use
Cutaway model of a backflow preventer
Mounted on demonstration panel; 500x400mm.

Technical Specification :

Technical Data:

Material: grey cast iron/steel
Dimensions and Weight
l x w x h: 500 x 370 x 400 mm
Weight: 15 kg.



Laboratory Instrument India