



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cutaway Model Changeover Valve

Product Code :
LBNY-0005-103000159



Description :

Functionality of moving parts preserved
Cutaway model of a changeover valve
Unit designed for desktop use
Mounted on demonstration panel; 500x400mm.

Technical Specification :

Technical Data:

Material: grey cast iron/steel
Flange fittings: DN 50
Dimensions And Weight
l x w x h: 500 x 370 x 400 mm
Weight: approx. 20 kg.



Laboratory Instrument India