



Email : sales@laboratoryinstrumentindia.com

Product Name :
Physics Lab Sets Optics

Product Code :
LBNY-0002-700010



Description :

Physics Lab Sets Optics

Technical Specification :

The Light box kit is specially designed for demonstrating basic concepts of reflection, refraction & colour mixing of different colours using various optical components.

The kit can be used by an individual student or it can be used for demonstration to a small group of students.

Physics Lab Sets Optics, Physics Lab Sets Optics Equipments, Physics Lab Sets Optics Tools, Physics Lab Sets Optics Tool Kits, Physics Lab Sets Optics Manufacturers, Physics Lab Sets Optics Suppliers from India, China, Kenya



Laboratory Instrument India