



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Assembly Exercise Shut Off Valve

**Product Code :**  
LBNY-0005-103000161



**Description :**

Shut-off valve consisting of housing, hand wheel, clamp cover, packing gland frame, taper and spindle  
Shut-off valve parts set, with associated set of small parts, in a sturdy case  
Spindle sealing based on the gland principle  
Assembly exercise for engineering training  
The kit forms part of the assembly, maintenance and repair practice line.

**Technical Specification :**

**Technical Data:**

Housing, hand wheel, clamp cover, packing gland  
Frame: grey cast iron  
Shut-off valve with flange connections:  
Stroke: 13mm  
Taper, seating ring, spindle, ring segment etc.:  
Stainless steel.



**Laboratory Instrument India**

