



Email : sales@laboratoryinstrumentindia.com

Product Name :
Assembly Exercise Shut Off Valve

Product Code :
LBNY-0005-103000161



Description :

Shut-off valve consisting of housing, hand wheel, clamp cover, packing gland frame, taper and spindle
Shut-off valve parts set, with associated set of small parts, in a sturdy case
Spindle sealing based on the gland principle
Assembly exercise for engineering training
The kit forms part of the assembly, maintenance and repair practice line.

Technical Specification :

Technical Data:

Housing, hand wheel, clamp cover, packing gland
Frame: grey cast iron
Shut-off valve with flange connections:
Stroke: 13mm
Taper, seating ring, spindle, ring segment etc.:
Stainless steel.



Laboratory Instrument India

