




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway Model Standard Venturi Meter

**Product Code :**  
LBNY-0005-103000163



**Description :**

Cutaway model of a standard venturi meter  
Unit designed for desktop use  
Mounted on demonstration panel; 500x400mm.

**Technical Specification :**

**Technical Data:**

Material: gun metal/steel  
Flange fittings: DN 50.



**Laboratory Instrument India**