



Email : sales@laboratoryinstrumentindia.com

Product Name :
Assembly Station Pipes and Valves and Fittings

Product Code :
LBNY-0005-103000167



Description :

Standard commercially available flanged fittings: shut-off valve, non-return valve, strainer, condensation Drain, inspection glass, ball valve, gate valve
Ball valve with cutting ring screw fitting
Assembly exercise for engineering training
Mobile frame for mounting of pipe network
Connection to water supply via hose with coupling
Piping network comprising pipe bends, elbows, T-pieces and transitions in nominal
Pipe connections via flanges or cutting ring screw fittings
Pressure vessel with manometer, connection via.

Technical Specification :

Technical Data:

Ball valve with cutting ring screw fitting
Brass, nickel-plated
Flange fittings
Grey cast iron
Nominal size: G1/2"
Nominal pressure: PN25
Manometer: 0...4bar.



Laboratory Instrument India