



Email : sales@laboratoryinstrumentindia.com

Product Name :
Maintenance Exercise Multi Stage Centrifugal Pump

Product Code :
LBNY-0005-103000168



Description :

Flow setting by manual stroke length adjustment (including during operation)
Manual drive with crank instead of a drive motor
Learning concept for maintenance and repair exercises on a single-diaphragm pump
Pump parts and tools housed in a tool box
Diaphragm and push rod directly linked.

Technical Specification :

Technical Data:
Diaphragm pump
Nominal stroke rate at 50Hz: 156min-1
Power consumption: max. 90W
Flow rate: 0...2,4L/h
Max. head: 100m
Pump materials
Pump body: Polypropylene (PP)
Double-ball valves: PP/glass fibre-reinforced plastic
Drive diaphragm: PTFE-lined
Valve balls: glassk
Valve seals: FPM.



Laboratory Instrument India