



Email : sales@laboratoryinstrumentindia.com

Product Name :
Physics Lab Sets Magnetism

Product Code :
LBNY-0002-700011



Description :

Physics Lab Sets Magnetism

Technical Specification :

Physics lab sets magnetism

This isotropic ferrite magnet kit is designed for kids and is the perfect training tool for science classrooms, helping teachers / professors who want to learn more about science / magnetism. Attracts attention to children. Experiment with these magnet kit sets to learn basic physics principles and explore physics and magnetism. Physics is a practical science. Good quality, proper physical activity and research not only motivates and enjoys but also stimulates students' observation and questions. Key to strengthening the integration of learning, clarification and theory.

A set of interesting magnets and accessories for physical science experiments.

Physics Lab Sets Magnetism, Physics Lab Sets Magnetism Equipments, Physics Lab Sets Magnetism Tools, Physics Lab Sets Magnetism Tool Kits, Physics Lab Sets Magnetism Manufacturers, Physics Lab Sets Magnetism Suppliers from India, China, Kenya



Laboratory Instrument India

