




Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Cutaway Model Straight Way Plug Valve

**Product Code :**  
LBNY-0005-103000171



**Description :**

Unit designed for desktop use  
Functionality of moving parts preserved  
Cutaway model of a straight-way plug valve  
Mounted on demonstration panel; 400x300mm.

**Technical Specification :**

**Technical Data:**

Material: gun metal/steel  
Fittings: 1 ".



**Laboratory Instrument India**