



Email : sales@laboratoryinstrumentindia.com

Product Name :
Maintenance Exercise In Line Centrifugal Pump

Product Code :
LBNY-0005-103000173



Description :

Pump shaft sealing with floating ring seal
Pump drive by 3-phase AC motor
Learning concept for maintenance and repair exercises on an in-line centrifugal pump
Pump parts and tools housed in a tool box
Enclosed pump impeller with 5 blades, designed for pure liquids.

Technical Specification :

Technical Data:
In-line centrifugal pump
Housing and impeller: grey cast iron
Max. flow rate: 19m³/h
Max. head: 16m
Power consumption: max. 750W
Speed: 2900min⁻¹
Drive motor
400V, 50Hz, 3 phases; or 230V, 60Hz, 3 phase.

Laboratory Instrument India

