



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Valves And Fittings Test Stand

Product Code :
LBNY-0005-103000174



Description :

Hand-operated piston pump to generate the test pressure, a return valve to relieve the system pressure, and a manometer for pressure measurement
2 different sizes of mounting flange with blank flange and flange seal
Test stand on which to mount industrial valves and fittings
Pressure testing of valves and fittings
Mobile frame with collector tray and ball valve to drain
Connection of pump and test flange via pressure hose
Test medium: water
Water storage tank.

Technical Specification :

Technical Data:
Piston pump with tank
Tank capacity: 12L
Manometer: 0...60bar
Test pressure: 0...60bar.

Laboratory Instrument India

