



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Circuit With Centrifugal Pump

Product Code :
LBNY-0005-103000185



Description :

Measuring objects: pump, flow meter, diaphragm valve
Pressure conditions at various measuring objects
Flow measurement using rotameter
Flow can be adjusted via valves
Centrifugal pump with 3 different speeds
Closed water circuit
Differential pressure measurement using electronic pressure meter
Annular chambers allow easy measurement of pressure.

Technical Specification :

Technical Data:

Tank
Volume: 8,5L
Pump:
Max. head: 6m
Max. flow rate: 5m/h
Max. power consumption: 70W
Three switching stages for speed selection
Measuring ranges
Differential pressure: -350...350mbar
Flow rate: 150...1600L/h.



Laboratory Instrument India