



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hydraulic Circuit With Centrifugal Pump

**Product Code :**  
LBNY-0005-103000185



**Description :**

Measuring objects: pump, flow meter, diaphragm valve  
Pressure conditions at various measuring objects  
Flow measurement using rotameter  
Flow can be adjusted via valves  
Centrifugal pump with 3 different speeds  
Closed water circuit  
Differential pressure measurement using electronic pressure meter  
Annular chambers allow easy measurement of pressure.

**Technical Specification :**

**Technical Data:**

Tank  
Volume: 8,5L  
Pump:  
Max. head: 6m  
Max. flow rate: 5m/h  
Max. power consumption: 70W  
Three switching stages for speed selection  
Measuring ranges  
Differential pressure: -350...350mbar  
Flow rate: 150...1600L/h.



**Laboratory Instrument India**