



Email : sales@laboratoryinstrumentindia.com

Product Name :
Base Unit For Turbines

Product Code :
LBNY-0005-103000190



Description :

Basic experiments on centrifugal pumps
Includes pump and transparent water tank
Together with the turbines: investigation of operating behaviour and Recording of turbine characteristics
Microprocessor-based measuring technique
Sensors for flow rate and pressure
Low air entry into circulating water ensured by damping plate inside the tank
Variable pump speed via frequency converter.

Technical Specification :

Technical Data:

Pump
Power consumption: 670W
Max. head: 35,4m
Max. flow rate: 70L/min
Water tank: approx. 15L
Measuring ranges
Pressure: -1...5bar
Flow rate: 3,9...50L/min.

Laboratory Instrument India

