



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hydraulic Ram Pumping Using Water Hammer

Product Code :
LBNY-0005-103000194



Description :

Formation and effect of water hammer
Waste valve with adjustable lift, closes cyclically due to flow force of the water
Tank with check valve and air volume is used as an air vessel
Pumping using water hammer
Elevated tank with variable pump head
Air volume in the air vessel is varied by vent valve
Fixed overflow tank is used as a water source, e.g. river, pool.

Technical Specification :

Technical Data:

Ram
Max. flow rate: 90L/h
Max. head 0,27m.



Laboratory Instrument India

