



Email : sales@laboratoryinstrumentindia.com

Product Name :
Wake Measurement

Product Code :
LBNY-0005-103000195



Description :

Measurement of drag via two-component force sensor
Removable spacer plate allows two wake rake positions for measurement
Investigation of the wake of a cylinder immersed in a flow
Accessory for the wind tunnel
Display of the pressures on the 16 tube manometers
Cylinder as drag body
Wake rake with 15 Pitot tubes detects total pressures.

Technical Specification :

Technical Data:

Cylinder
Height: 100mm
Inner diameter: 20mm
Height with supporting rod: 290mm
Wake rake
15 Pitot tubes
Distance between the Pitot tubes: 3mm
Inner diameter: 1,1mm
Outer diameter: 2mm.



Laboratory Instrument India