



Email : sales@laboratoryinstrumentindia.com

Product Name :
Aerofoil Spring Mounted

Product Code :
LBNY-0005-103000205



Description :

Profile: Symmetrical
Aerofoil for examination of self-starting vibration

Technical Specification :

Technical Data:

Aerofoil
Mass moment of inertia: $1,07 \cdot 10^{-4}$ kgm
L x W x H 200x100x15mm
Weight: 0,157kg
Spring rigidity
Torsional rigidity: 0,07...0,28 Nm/rad
Total rigidity, transverse: 216 N/m
Single spring: 27 N/m.



Laboratory Instrument India

