




Email : sales@laboratoryinstrumentindia.com

Product Name :
Visualisation Of Streamlines

Product Code :
LBNY-0005-103000206



Description :

Visualisation Of streamlines by using fog
Vertical measuring section with transparent front plate and black background
Aerofoil with adjustable angle of attack
Scale for displaying the angle of attack
Fog generator, operation with non-toxic and water-soluble fog fluid
Three models for insertion into the wind tunnel
Accessories for Aerodynamics Trainer.

Technical Specification :

Technical Data:

Models
Aerofoil, adjustable angle of attack
Cylinder
Orifice plate
Measuring Section
Cross-section in the viewing area: 252x42mm.

Laboratory Instrument India

