



Email : sales@laboratoryinstrumentindia.com

Product Name :
Cavitation In Pumps

Product Code :
LBNY-0005-103000208



Description :

Open impeller to observe the blades during operation
Continuously adjustable pump speed via frequency converter
Visualisation of cavitation in centrifugal pumps
Transparent pump housing and pipe at the inlet side
Display of the pressures at inlet and outlet side of the pump via manometers
Closed water circuit with tank and temperature display
Digital display of speed, water temperature in the return and flow rate
Temperature control via heater and external cooling via water supply
Flow measurement using rotameter.

Technical Specification :

Technical Data:

Centrifugal pump with drive motor
Tank: 20L
Power consumption: 0,37kW
Max. head: 13m
Max. flow rate: 70L/min
Speed: 500...3300min-1.

Laboratory Instrument India

