



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Networks

Product Code :
LBNY-0005-103000212



Description :

Investigation of different pipe networks
Construction of pipe networks from pipe sections and various piping elements
Calibration of pipe sections
Parallel and series connection of pipe sections
Differential pressure measurement with twin tube manometers and differential pressure manometer
Five pre-installed pipe sections with different diameters
Panel for piping elements
Construction of a closed circular pipeline
Flow rate measurement with measuring tank (can be shut off), stopwatch and level indicator.

Technical Specification :

Technical Data:

Pump
Max. head: 7,6m
Max. flow rate: 9m/h
Power consumption: 250W
Pipe network, max. flow rate: 4,8m/h
Pipe sections, length 700mm each
1x: 25x1,9mm
2x: 20x1,5mm
2x: 16x1,2mm
Twin tube manometers

Tank for water: 180L
Tank for flow rate measurement
Small measuring range: 10L
Large measuring range: 60L
Stopwatch: 1/100s.



Laboratory Instrument India