




Email : sales@laboratoryinstrumentindia.com

Product Name :
Pump And Valves And Fittings Test Stand

Product Code :
LBNY-0005-103000215



Description :

Trainer for testing various valves and fittings
Pressure measuring points upstream and downstream of test device for differential pressure Manometer with pressure switch
Installation of the test fitting in a pipe section of variable length
Centrifugal pump with variable speed via frequency converter
Digital displays for flow rate, power output, speed, position of control valve
Fine pressure regulator adjusts compressed air pressure
Manometers at inlet and outlet of centrifugal pump
Tank cover as collecting tray under test device.

Technical Specification :

Technical Data:

Centrifugal pump
Speed: 1450...2900min⁻¹
Max. flow rate: 72m³/h
Max. head: 26,5m
Power consumption: 4kW
Supply tank with cover: capacity: 400L
Test valves and fittings:
Shut-off valve
Ball valve with pneumatic drive
Safety valve 1", 1,5bar

Wedge gate valve
Electric control valve
Butterfly valve.



Laboratory Instrument India