



Email : sales@laboratoryinstrumentindia.com

Product Name :

Orifice Plate Flow Meter With Transducer

Product Code :

LBNY-0005-103000216



Description :

Orifice Plate Flow Meter with electronic differential pressure transducer for flow rate measurement as accessory
Vertical and horizontal installation possible
Connections to facilitate pressure loss measurement
Connections to supply auxiliary power via
Operation based on the differential pressure method
Display indicating differential pressure.

Technical Specification :

Technical Data:

Orifice plate
Material: brass
Transducer measuring range: 0...500mbar
Diameter: 18,5mm, with 45 chamfer
Pipe connections: DN 32
Auxiliary power: 24VDC.



Laboratory Instrument India

