



Email : sales@laboratoryinstrumentindia.com

Product Name :
Heat Transfer At A Cylinder In Transverse Flow

Product Code :
LBNY-0005-103000220



Description :

Copper cylinder with integrated temperature sensor
Convective heat transfer of a cylinder in an air-flow tube
Oven keeps the temperature constant
Cylinder is heated in the oven to approximately 120°C
Display and analysis of the measured values using the software.

Technical Specification :

Technical Data:
Oven
Heating power: 100W
Cylinder
Outer diameter: 20mm
Material: copper
Length: 120mm
Measuring range
Temperature: 0...200°C.

Laboratory Instrument India

