




Email : sales@laboratoryinstrumentindia.com

Product Name :
Electronic Total Pressure Sensor

Product Code :
LBNY-0005-103000221



Description :

Vertical adjustment of the Pitot tube
Electronic detection of the position
Electronic total pressure sensor for measurement of static and dynamic pressure
Display and analysis of the measured values using the software
In conjunction with measurement of the total pressure in the wake of a cylinder.

Technical Specification :

Technical Data:

Pitot tube
Inner diameter: 0,41mm
Outer diameter: 0,71mm
Vertical adjustment: 0...130mm
Dimensions and Weight
LxWxH (retracted): 120x75x350mm
Weight: approx. 1kg.



Laboratory Instrument India

