




Email : sales@laboratoryinstrumentindia.com

Product Name :
Measurement Of Boundary Layers

Product Code :
LBNY-0005-103000225



Description :

Vertically sliding Pitot tube can be set to 3 positions along the plate in the measuring section
Display of static and total pressure on the tube manometer
Additional measuring point for measuring the static pressure
Investigation of the boundary layer on a flat plate
Accessory for experimental plant
Plate leading edge with chamfer
Pitot tube for measuring the total pressure.

Technical Specification :

Technical Data:

Flat plate
LxWxH: 420x80x8mm
15 chamfer facing the incident flow
Pitot tube
Inner diameter: 0,6mm
Vertically sliding: 0...18mm
Measuring section with 3 positions along the
Plate: 10mm, 210mm and 410mm from the leading.

Laboratory Instrument India

