



Email : sales@laboratoryinstrumentindia.com

Product Name :
Free Jet

Product Code :
LBNY-0005-103000226



Description :

Determination of the pressure loss in the flow outlet into stationary surroundings
16 tube manometers for displaying the pressures is recommended
Pitot tube, can be moved horizontally and vertically
Accessory for aerodynamics trainer.

Technical Specification :

Technical Data:

Pitot tube, sliding
Inner diameter: 2mm
Vertical: 0...700mm
Horizontal: -150...150mm.
Pipe socket, plastic
Inner diameter: 54mm.



Laboratory Instrument India

