



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Free Jet

**Product Code :**  
LBNY-0005-103000226



**Description :**

Determination of the pressure loss in the flow outlet into stationary surroundings  
16 tube manometers for displaying the pressures is recommended  
Pitot tube, can be moved horizontally and vertically  
Accessory for aerodynamics trainer.

**Technical Specification :**

**Technical Data:**

Pitot tube, sliding  
Inner diameter: 2mm  
Vertical: 0...700mm  
Horizontal: -150...150mm.  
Pipe socket, plastic  
Inner diameter: 54mm.



**Laboratory Instrument India**

