



Email : sales@laboratoryinstrumentindia.com

Product Name :
Paddle Wheel Flow Meter

Product Code :
LBNY-0005-103000228



Description :

Non-contact recording of rotation speed by ring magnet and Hall sensor
Display indicating flow rate
Paddle wheel flow meter as accessory for trainer
Vertical and horizontal installation possible
Connections to supply auxiliary power via
Connections to facilitate pressure loss measurement.

Technical Specification :

Technical Data:
Auxiliary power: 24VDC
Measuring range: 2...48L/min.



Laboratory Instrument India