



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Paddle Wheel Flow Meter

**Product Code :**  
LBNY-0005-103000228



**Description :**

Non-contact recording of rotation speed by ring magnet and Hall sensor  
Display indicating flow rate  
Paddle wheel flow meter as accessory for trainer  
Vertical and horizontal installation possible  
Connections to supply auxiliary power via  
Connections to facilitate pressure loss measurement.

**Technical Specification :**

**Technical Data:**

Auxiliary power: 24VDC  
Measuring range: 2...48L/min.



**Laboratory Instrument India**