



Email : sales@laboratoryinstrumentindia.com

Product Name :
Friction Losses In Pipe Elements

Product Code :
LBNY-0005-103000231



Description :

Measuring points on straight pipe sections for pressure measurement
Display and analysis of the measured values using
Investigate friction losses in various pipe elements
Measuring points at pipe inlets for velocity measurement
Two interchangeable pipe inlets: sudden enlargement or uninterrupted air inlet
Combine different pipes from pipe elements.

Technical Specification :

Technical Data:
Straight pipe section with flange
Inner diameter: $d=53,6\text{mm}$
Length: 1235mm
Straight pipe section
Length: 991mm
Pipe inlets
Sudden enlargement: $d=35\dots53,6\text{mm}$
Uninterrupted air inlet: radius= 22mm
90 pipe bend
Inner diameter: $d=53,6\text{mm}$
90 pipe angle
Radius: $2xd$
Inner diameter: $d=53,6\text{mm}$.



Laboratory Instrument India