



Email : sales@laboratoryinstrumentindia.com

Product Name :
Methods Of Flow Measurement

Product Code :
LBNY-0005-103000234



Description :

Different methods of flow rate measurement
Measurement of the total pressure with Pitot tube
6 tube manometers to determine the pressure distribution in Venturi nozzle, orifice plate flow meter and measuring nozzle
Measuring instruments: orifice plate flow meter/measuring nozzle, Venturi nozzle and rotameter.

Technical Specification :

Technical Data:
Venturi nozzle: A=84...338mm
Angle at the outlet: 4
Angle at the inlet: 10,5
Measuring nozzle: diameter=18,5mm
Orifice plate flow meter: diameter=14mm
6 tube manometers: 390mmWC
Rotameter: max. 1700L/h.



Laboratory Instrument India

