



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Methods Of Flow Measurement

**Product Code :**  
LBNY-0005-103000234



**Description :**

Different methods of flow rate measurement  
Measurement of the total pressure with Pitot tube  
6 tube manometers to determine the pressure distribution in Venturi nozzle, orifice plate flow meter and measuring nozzle  
Measuring instruments: orifice plate flow meter/measuring nozzle, Venturi nozzle and rotameter.

**Technical Specification :**

**Technical Data:**  
Venturi nozzle: A=84...338mm  
Angle at the outlet: 4  
Angle at the inlet: 10,5  
Measuring nozzle: diameter=18,5mm  
Orifice plate flow meter: diameter=14mm  
6 tube manometers: 390mmWC  
Rotameter: max. 1700L/h.



**Laboratory Instrument India**

