



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Series And Parallel Configuration Pumps

**Product Code :**  
LBNY-0005-103000235



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Investigation of series and parallel configuration of pumps  
Two identical centrifugal pumps  
Transparent tank as intake tank  
Overflow in the tank ensures constant suction head  
Ball valves used to switch between series and parallel operation  
Manometers at inlet and outlet of each pump.

**Technical Specification :**

**Technical Data:**  
2x centrifugal pump  
Power consumption: 370W  
Max. flow rate: 21L/min  
Max. head: 12m  
Tank: 13L  
Pipes and pipe connections: PVC  
Measuring ranges  
Pressure (inlet): 2x -1...1,5bar  
Pressure (outlet): 3x 0...2,5bar

---

# Laboratory Instrument India

