




Email : sales@laboratoryinstrumentindia.com

Product Name :
Vertical Flow From A Tank

Product Code :
LBNY-0005-103000239



Description :

Measuring device for determining the jet diameter
Determining the contraction coefficient for different contours and diameters
Pitot tube for determining the total pressure
Pressure display on twin tube manometers
Tank with adjustable overflow
5 interchangeable inserts with different contours
Study of pressure losses in vertical flows from tanks.

Technical Specification :

Technical Data:

Tank
Capacity: approx. 13L
Max. flow rate: 14L/min
Overflow height: max. 400mm.



Laboratory Instrument India

