




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Vertical Flow From A Tank

**Product Code :**  
LBNY-0005-103000239



**Description :**

Measuring device for determining the jet diameter  
Determining the contraction coefficient for different contours and diameters  
Pitot tube for determining the total pressure  
Pressure display on twin tube manometers  
Tank with adjustable overflow  
5 interchangeable inserts with different contours  
Study of pressure losses in vertical flows from tanks.

**Technical Specification :**

**Technical Data:**

Tank  
Capacity: approx. 13L  
Max. flow rate: 14L/min  
Overflow height: max. 400mm.



**Laboratory Instrument India**

