



Email : sales@laboratoryinstrumentindia.com

Product Name :
Flow Visualisation Apparatus

Product Code :
LBNY-0005-103000240



Description :

Flow velocity, water supply and water drain in sources/sinks as well as dosage of the contrast medium Can be adjusted by using valves
Visualisation of streamlines
Water as flowing medium and ink as contrast medium
Bottom plate with water connections for generating sources/sinks
Sources/sinks can be combined as required
Rubber plate for creating your own models included
Upper glass plate, hinged for interchanging models
Different drag bodies and changes in cross-section included.

Technical Specification :

Technical Data:

Flow chamber contains two plates
Upper plate made of glass
Distance between the plates: 2mm
10 drag bodies and changes in cross-section
Bottom glass plate with four water connections for Sources/sinks
Size experiment area: LxW: 400x280mm
Rubber plate for your own models
Thickness: 2mm
LxH: 300x400mm

Tank for contrast medium: 500mL
Injection of the contrast medium
15 holes.



Laboratory Instrument India