



Email : sales@laboratoryinstrumentindia.com

Product Name :
Pipe Friction For Laminar And Turbulent Flow

Product Code :
LBNY-0005-103000241



Description :

Transparent tank with overflow ensures constant water inlet pressure in the pipe section for experiments with laminar flow
Investigation of the pipe friction in laminar or turbulent flow
Flow rate adjustment via valves
Twin tube manometers for measurements in laminar flow
Dial-gauge manometer for measurements in turbulent flow
Water supply via or via laboratory supply for experiments with turbulent flow.

Technical Specification :

Technical Data:

Tank: 5L
Pipe Section
Inside diameter: 3mm
Length: 400mm
Measuring ranges
Differential pressure (dial-gauge manometer): 0...0,4bar
Differential pressure (twin tube manometers): 2x 370mmWC.

Laboratory Instrument India

