



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Pipe Friction For Laminar And Turbulent Flow

**Product Code :**  
LBNY-0005-103000241



**Description :**

Transparent tank with overflow ensures constant water inlet pressure in the pipe section for experiments with laminar flow  
Investigation of the pipe friction in laminar or turbulent flow  
Flow rate adjustment via valves  
Twin tube manometers for measurements in laminar flow  
Dial-gauge manometer for measurements in turbulent flow  
Water supply via or via laboratory supply for experiments with turbulent flow.

**Technical Specification :**

**Technical Data:**

Tank: 5L  
Pipe Section  
Inside diameter: 3mm  
Length: 400mm  
Measuring ranges  
Differential pressure (dial-gauge manometer): 0...0,4bar  
Differential pressure (twin tube manometers): 2x 370mmWC.

---

# Laboratory Instrument India

